SAFE-ST NANO SPITTOON FOAM CLEANER & DISINFECTANT

- Nano Foam Cleanser for disinfecting spittoons
- Ready-to-use nano foam cleanser for cleaning and disinfecting spittoons

Photo

CHARACTERISTICS:

- Active cleaning in aerosol form
- Aldehyde-free
- Alcohol based
- Ready to use
- Easy to apply foam (no contact)

PRODUCT DESCRIPTION:

The cleanser is active against a broad spectrum of bacteria and yeasts.

The spittoon disinfectant is a mild and ready-to-use bactericidal special foam. It is suitable for disinfection, cleaning and maintenance of spittoons in dental treatment units and reliably removes lime scale deposits, blood and residue from plaque disclosing tablets with the help of an active cleaning system. Do not apply to alcohol-sensitive surfaces.

APPLICATION:

The spittoon disinfectant is used undiluted. Briefly shake the can before use. To open the safety closure, press on the smooth surfaces and twist. Cover the entire surface of the spittoon with about 20 ml of the cleaner and leave the solution to work. If needed, spread the cleaning solution around with a sponge or brush. After at least one minute application time, rinse off with water. During use, hold the can with the spray nozzle facing downwards and activate the valve. Spray the foam all around the spittoon and leave the solution to work. After at least one minute application, either mechan-

ically clean the spittoon or rinse it with water. Please note that the foam will collapse within the one minute of application.

It is recommended to use the foam after each contamination.

APPLICATION TIME:

Bactericidal effect: Bactericide as defined by EN 13697 (Staphylococcus aureus, Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa), low load: undiluted for 1 minute

Yeasticidal effect: Candida albicans as defined in EN 13697 (low load): undiluted for 1 minute

Tuberculocidal effect: Tuberculosis inactivation (M. terreae) as defined in EN 14348 (low load): undiluted for 5 minutes

COMPOSITION:

Bis-(3-aminopropyl)-dodecylamine, ethanol, excipient, propellant, water. pH: approx. 10.5

SAFE-ST NANO SPITTOON FOAM CLEANER & DISINFECTANT

- Nano Foam Cleanser for disinfecting spittoons
- Ready-to-use nano foam cleanser for cleaning and disinfecting spittoons

Photo

SPECIAL INFORMATION:

Safety Precautions: Pressurized container. Keep out of the reach of children. Extremely flammable aerosol. Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Do not spray on an open flame or other ignition source. Pressurized container: Do not pierce or burn, even after use. Avoid breathing dust/fume/gas/mist/vapours/spray. Protect from sunlight. Do not expose to temperatures exceeding 50°C/ 122°F. Pressurized container: may burst if heated. Contains: propane, butane, ethanol, N- (3-aminopropyl) -N-1,3-diamine.

H226 Flammable liquid and vapour. H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. H412 Harmful to aquatic life with long lasting effects. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations. 648/2004/EC: < 5% nonionic surfactants, < 5% EDTA.

PRODUCT STATUS/LISTINGS:

Conform to the EU guideline 93/42/EC for medical devices.

AVAILABLE DELIVERY FORMS:

Itemnumber	Packaging
9797742	400ml Aluminium can (aerosol dispenser for foam application)

Use disinfectants safely. Always read the label and product information before use.

SAFE-ST NANO SPITTOON FOAM CLEANER & DISINFECTANT

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