



Light- Activated Temporary Material

SECTION I - Product/ Company Identification

Product Trade Name: Tempit® L/C Light- Activated Temporary Material for Inlays & Onlays

Part (Item) Numbers: 310063
Manufacturer: CENTRIX INC
Address: 770 RIVER ROAD
City, State, Zip: SHELTON, CT 06484

Emergency Telephone Number: (800) 235-5862 USA & CANADA Information Telephone Number: (800) 235-5862 (203) 929-5582

SECTION II - Hazardous Ingredients / Identity Information

Hazardous Components:	%	CAS NO.	SYMBOL	R	OTHER LIMITS
HEMA	6-10%	868-77-9	Xi	36-38	
UDMA	10-14%	41137-60-4	Xi	36-38	
Carbomer 941	40-45%	US-Pharmacop	Xi	36-38	

SECTION III - Hazards Identification

Designation of risk: Xi, irritant Special details of risk (R-categories): 36/37/38

Detrimental effects/symptoms: Irritant in case of eye, skin and inhalation contact.

SECTION IV - First Aid Measures

Skin: Rinse thoroughly with plenty of water. Ingestion: Contact immediately a physician.

Inhalation: Fresh air, in case of indisposition or nausea contact a doctor.

Eyes: Rinse thoroughly with water, contact a doctor.

Information for the attending physician: Material for dental or laboratory use only.

SECTION V - Fire-Fighting Measures

Suitable extinguishing media: Foam, carbon dioxide, powder fire extinguishers.

Forbidden extinguishing media: Water. Unusual fire and explosion hazards: None.

Special protection equipment: No data available.

SECTION VI - Accidental Release Measures

Personal precautions: Avoid eye and/or skin contact.

Environmental precautions: Do not let the material get in contact with sewage systems.

Reclaiming methods: Take up mechanically, use comman absorbents.

SECTION VII - Handling and Storage

Handling: Only by qualified staff in dental offices or dental laboratories.

Storage: Tightly closed at room temperature.
Containers: Store in good ventilated roomsand areas.

Incompatibilities: Do not store together with food. Other precautions: Keep out of the reach of children.

SECTION VIII - Exposure Controls / Personal Protection

General measures: None.

Limits:

Respiratory protection: Do not breeth or inhale vapour or dust.

Eye protection: Goggles/Eyeshields.
Skin protection: Overalls, rubber gloves.

Hands protection: Overalls, rubber gloves. And wash hands after use.

Other equipment:

SECTION IX - Physical and Chemical Properties

Appearance: Yellow to White Paste

Explosion limits: N/A

Odor: Methacrylate-like

Water solubility: mg/1 Boiling point: N/A

Evaporation rate:

Flash point: N/A

SECTION X - Stability and Reactivity

Stability:

Conditions to avoid: Light exposure, heat.

Materials to avoid: Radical initiators, e.g., peroxides Hazardous decomposition products: HCN, cyanides

Further particulars: Material cures in case of light exposition.

SECTION XI - Toxicological Information

SECTION XII - Ecological Information

SECTION XIII - Disposal Considerations

In accordance with local authorities.

SECTION XIV - Transport Information

SECTION XV - Regulatory Information

Hazard Classification Symbol:

Ingredients determining the label: The information is this MSDS does not constitute an assessment of workplace risk.

SECTION XVI - Other Information

Handling of product exclusively by dental specialists according to instructions for use.

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY.

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.



Making Dentistry Easier.[™]

CENTRIX INCORPORATED

770 River Road, Shelton, CT 06484 USA (800) 235-5862 Tel: (203) 929-5582 Fax: (203) 929-6804



0510

CENTRIX EUROPE

Hauptstrasse 64 B, 65719 Hofheim, Germany

©Centrix, Inc. 2008 Date Revised 10/19/2010